

How to choose the “right” MINI PAK’R™ cushion for your packing job? Ask Art the Engineer!



*Every cushion
has a purpose.
Ask me!*

Use NOVUS Quilt-Large Cushions if you are:



Shipping large dishware, candles or other fragile items.
Wrap and cushion with Quilt Large.

Use NOVUS Quilt-Small Cushions to:



Wrapping decorative objects, glasses,
other breakable items.



Interleaving between stacks of plates
or picture frames.

Use NOVUS Double Cushions to:



Top fill a carton quickly with fewer cushions. Fill voids
and block and brace items in a carton.

Still not sure which is “right” MINIPAK’R™ cushion for your packing job?



“Ask Art the Engineer!”

Use NOVUS Supertube Cushions for:



Shipping cosmetics, fragile glassware, electronics. Use to line the carton, to wrap corners of boxes or books, and to wrap breakable items.

Use CELL-O 8” x 8” Cushions to:



Fill large voids, for top fill, and to block and braces boxes or books in a carton.

Use CELL-O 4” x 8” Cushions for:



Blocking and bracing items in a carton. Also to fill small voids and for top fill.

www.fpintl.com
877-594-5160

FP
INTERNATIONAL®
INVENTING THE FUTURE OF PACKAGING